## DL 2000-G DNA Ladder

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产品名称	规格	Cat. No.
DL 2000-G DNA Ladder	250 $\mu$ l	MD1006-G
DL 2000-G DNA Ladder	1.25 ml	MD1106-G

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## STORAGE

Store at room temperature (15 - 25°C). For long term storage, store at -20°C.


## TECHNICAL SUPPORT

TEL: 400-0099-857

E-MAIL: [technical@simgen.cn](mailto:technical@simgen.cn)

## INTRODUCTION

The DL 2000-G DNA Ladder consists of 6 individual DNA fragments ranging from 100 bp to 2000 bp and genomic DNA. DL 2000-G Ladder DNA Marker contains 1  $\times$  Loading Buffer, users can apply 5 - 10  $\mu$ l Marker in agarose gel electrophoresis directly.

## PRECAUTION

1. Clear bands can be obtained by applying 5  $\mu$ l Ladder when the lane width is less than 5 mm. If the lane is wider, loading volume of Ladder should be increased appropriately.
2. For DNA electrophoresis, agarose purity is of great significance to DNA band definition. Therefore, agarose with good quality should be used and gel concentration of 2 - 2.5 % is recommended.
3. During agarose electrophoresis, the concentration of agarose is closely associated with the separation of DNA fragments. High agarose concentration is ideal for the separation of the short DNA fragments. While low agarose concentration is ideal to separate the long DNA fragments.
4. The red and yellow tracking dye in DL 2000-G DNA Ladder will not weakened the DNA bands under UV light, better than bromophenol blue and xylene cyanol FF.