DL8000 Ladder

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| 产品名称 | 产品规格 | Cat. No. |
|---------------|------------------------|----------|
| DL8000 Ladder | 250 μ l | MD1015 |
| DL8000 Ladder | 250 μ l \times 5 | MD1115 |

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STORAGE

The product can be stored at room temperature (0-30 $^{\circ}$ C) for 3 years. If it is not used for a long time, please store it at - 20 $^{\circ}$ C to prevent moisture evaporation, and the validity period can be extended.

TECHNICAL SUPPORT


TEL: 400-0099-857

E-MAIL: technical@simgen.cn

INTRODUCTION

DL8000 Ladder is composed of 6 individual DNA fragments, presenting 8k, 5k, 3k, 1.5k, 1k, 500 bp sharp bands respectively. DL8000 Ladder contains 1 \times Loading Buffer, users can apply 5 - 10 μ l in agarose gel electrophoresis directly.

The red and yellow tracking dye in DL8000 Ladder will not weakened the DNA bands under UV light, better than bromophenol blue and xylene cyanol FF.

PRECAUTION

1. Clear bands can be obtained by applying 5 μ l DNA Ladder when the lane width is less than 5 mm. If the lane is wider, loading volume of DNA Ladder should be increased appropriately.
2. For DNA electrophoresis, agarose purity is of great significance to DNA band definition. Therefore, agarose with good quality should be used and gel concentration of 1% is recommended.
3. During agarose electrophoresis, the concentration of agarose is closely associated with the separation of DNA fragments. High agarose concentration is ideal for the separation of the short DNA fragments. While low agarose concentration is ideal to separate the long DNA fragments.