Supercoiled DNA Ladder

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产品名称	规格	Cat. No.
Supercoiled DNA Ladder	250 μ l	MD1014

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Supercoiled DNA Ladder

MATERIALS PROVIDED

Components	Specification	Cat. No.
Supercoiled DNA Ladder	250 µl	MD1014

STORAGE

Store at room temperature (15 - 25°C). If the product is not used for a long period of time, please store at -20°C to prevent the evaporation of water.


TECHNICAL SUPPORT

TEL: 400-0099-857

E-MAIL: technical@simgen.cn

INTRODUCTION

Supercoiled DNA Ladder is composed of 8 individual supercoiled plasmid DNA, presenting 15k, 10k, 8k, 6k, 5k, 4k, 3k, 2k bp sharp bands respectively. Supercoiled DNA Ladder contains 1×Loading Buffer, users can apply 5 - 10 µl Marker in agarose gel electrophoresis directly.

PRECAUTION

1. Clear bands can be obtained by applying 5 µl Ladder when the lane width is less than 5 mm. If the lane is wider, loading volume of Ladder should be increased appropriately.
2. For DNA electrophoresis, agarose purity is of great significance to DNA band definition. Therefore, agarose with good quality should be used and gel concentration of 0.7 - 1 % is recommended.
3. During agarose electrophoresis, the concentration of agarose is closely associated with the separation of DNA fragments. High agarose concentration is ideal for the separation of the short DNA fragments. While low agarose concentration is ideal to separate the long DNA fragments.
4. The red and yellow tracking dye in Supercoiled DNA Ladder will not weakened the DNA bands under UV light, better than bromophenol blue and xylene cyanol FF.