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3. Add 350 µl Buffer N8, mix gently and thoroughly by inverting the tube until the color of precipitate changes from blue to light yellow completely.

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* Continue to invert gently 4 - 6 times after the color of precipitate changes from blue to light yellow completely if the bacteria culture medium used was more than 3 ml.

4. Centrifuge at full speed ($\geq 12,000$ rpm) for 2 min.

5. Place a spin column in a 2 ml collection tube (provided), transfer the supernatant into the spin column, close the lid and centrifuge at 12,000 rpm for 30 s.

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6. Discard the filtrate. Place the spin column back to the collection tube. Add 800 µl Buffer W2 to the spin column. Closed the lid and centrifuge at 12,000 rpm for 30 s.

* Repeat step 6 is recommend to ensure high quality of plasmid, if plasmid will be used for cell transfection.

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* For endA+ host bacteria (e.g. wild-type E.coli, HB101, BL21), washing the spin column by adding Buffer W1 (Order separately, Cat. No. B025-01) is necessary to remove trace nuclease activity. Buffer W1 is not provided, can be purchased from Simgen separately.

7. Discard the filtrate. Place the spin column back to the collection tube. Centrifuge at full speed ($\geq 12,000$ rpm) for 1 min.

* Do not omit this step, as this will cause problems in downstream applications due to the residual ethanol in the eluate.

8. Discard the collection tube and place the spin column in a new 1.5 ml microcentrifuge tube. To elute DNA, add 60~100 µl Buffer E to the center of the column, close the lid, let stand for 1 min, and centrifuge at 12,000 rpm for 30 s.

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9. Discard the spin column. Eluted DNA can be used for a variety of molecular biology experiments or store at -20°C for later use.

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Protocol: Plasmid DNA Purification Using the Plasmid DNA Mini Kit

1. Harvest 1 - 5 ml overnight cultures of bacteria cells in LB (Luria-Bertani) medium by centrifuging at 12,000 rpm for 30 s, discard the medium. Fully suspend pelleted bacterial cells in 250 µl Buffer I with addition of RNase A.

* The pellet could be suspended completely by vortexing or pipetting up and down until no cell clumps left. Otherwise it would seriously affect the final plasmid DNA yield.

2. Add 250 µl Buffer II and mix thoroughly by gently inverting the tube 4 - 6 times.

* Before using Buffer II, make sure that no visible precipitate exists in the solution. After applying Buffer II, tighten the cap to avoid prolonged contact with air.

* Do not allow the lysis reaction to proceed for more than 5 minutes.

* Do not vortex, as this will result in shearing of genomic DNA.

* The solution should be viscous and slightly clear when the bacterial is lysed sufficiently; if too much bacteria was used, continue inverting the tube until the solution becomes viscous and slightly clear.

* Do not allow this step to proceed for more than 5 minutes.

3. Add 350 µl Buffer III, mix gently and thoroughly by inverting the tube until the color of precipitate changes from blue to light yellow completely.

* Do not vortex, as this will result in shearing of genomic DNA.

4. Centrifuge at full speed ($\geq 12,000$ rpm) for 10 min.

5. Place a spin column in a 2 ml collection tube (provided), transfer the supernatant in step 4 into the spin column, close the lid and centrifuge at 12,000 rpm for 30 s.

* Transfer the supernatant into the spin column immediately after centrifugation, stay for a long time may cause the precipitate to float, which makes it hard to operate in following procedure.

6. Discard the filtrate. Place the spin column back to the collection tube. Add 500 µl Buffer W1 to the spin column. Close the lid and centrifuge at 12,000 rpm for 30 s.

* This step is necessary to remove trace nuclease activity when using endA+ strains such as HB101, BL21. Host strains such as DH5 α do not require this additional wash step.

7. Discard the filtrate. Place the spin column back to the collection tube. Add 700 µl Buffer W2 to the spin column. Closed the lid and centrifuge at 12,000 rpm for 30 s.

* Ensure ethanol has been added into Buffer W2.

8. Discard the filtrate. Place the spin column back to the collection tube. Centrifuge at full speed ($\geq 12,000$ rpm) for 1 min.

* Do not omit this step, otherwise, it may cause problems in downstream applications due to the residual ethanol in the eluate.

9. Discard the collection tube and place the spin column in a new 1.5 ml microcentrifuge tube. To elute DNA, add 60~100 µl Buffer E to the center of each column, close the lid, let stand for 1 min, and centrifuge at 12,000 rpm for 30 s.

* The deionized water can also be used to elute the DNA, but the pH should be ensured at 7.0 - 8.5, otherwise it may affect the elution efficiency.

10. Discard the spin column. Eluted DNA can be used for a variety of molecular biology experiments or store at -20°C for later use.

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Protocol: Plasmid DNA Purification Using the Plasmid DNA Double Mini Kit

1. Harvest 5–15 ml overnight cultures of bacteria cells in LB (Luria-Bertani) in a 2 ml microcentrifuge tube (not provided), discard the medium. Fully suspend pelleted bacterial cells in 500 µl Buffer I with addition of RNase A.

* Recommended: Collect more than 10 ml of bacterial(except rich culture medium such as 2×YT) to increase the recovery of the final plasmid DNA. Bacteria can be collected in 15 ml or 50 ml centrifuge tubes, then suspended bacteria in Buffer I and transferred to 2 ml microcentrifuge tube.

* Please use a 2 ml microcentrifuge tube. The volume of the 1.5 ml microcentrifuge tube is too small to perform the follow steps.

* The pellet could be suspended completely by vortexing or pipetting up and down until no cell clumps left. Otherwise it would seriously affect the final plasmid DNA yield.

2. Add 500 µl Buffer II and mix thoroughly by gently inverting the tube 4 - 6 times.

* Before using Buffer II, make sure that no visible precipitate exists in the solution. After applying Buffer II, tighten the cap to avoid prolonged contact with air.

* Do not allow the lysis reaction to proceed for more than 5 minutes.

* Do not vortex, as this will result in shearing of genomic DNA.

* The solution should be viscous and slightly clear when the bacterial is lysed sufficiently; if too much bacteria was used, continue inverting the tube until the solution becomes viscous and slightly clear.

* Do not allow this step to proceed for more than 5 minutes.

3. Add 700 µl Buffer III, mix gently and thoroughly by inverting the tube until the color of precipitate changes from blue to light yellow completely.

* Do not vortex, as this will result in shearing of genomic DNA.

4. Centrifuge at full speed ($\geq 12,000$ rpm) for 10 min.

5. Place a spin column in a 2 ml collection tube (provided), transfer 800 µl supernatant from step 4 into the spin column, close the lid and centrifuge at 12,000 rpm for 30 s.

* If precipitate was suspended in the supernatant, repeat step 4 once again.

6. Discard the filtrate. Place the spin column back to the collection tube, and transfer remaining supernatant into the spin column, close the lid and centrifuge at 12,000 rpm for 30 s.

7. Discard the filtrate. Place the spin column back to the collection tube. Add 500 µl Buffer W1 to the spin column. Close the lid and centrifuge at 12,000 rpm for 30 s.

8. Discard the filtrate. Place the spin column back to the collection tube. Add 700 µl Buffer W2 to the spin column. Close the lid and centrifuge at 12,000 rpm for 30 s.

* Ensure ethanol has been added into Buffer W2.

9. Repeat step 8 once.

10. Discard the filtrate. Place the spin column back to the collection tube. Centrifuge at full speed ($\geq 12,000$ rpm) for 1 min.

* Do not omit this step, otherwise, it may cause problems in downstream applications due to the residual ethanol in the eluate.

11. Discard the collection tube and place the spin column in a new 1.5 ml microcentrifuge tube. To elute DNA, add 60–200 µl Buffer E to the center of each column, close the lid, let stand for 1 min, and centrifuge at 12,000 rpm for 30 s.

* The deionized water can also be used to elute the DNA, but the pH should be ensured at 7.0 - 8.5, otherwise it may affect the elution efficiency.

12. Discard the spin column. Eluted DNA can be used for a variety of molecular biology experiments immediately or store at -20°C for later use.

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1. Zhang G, Wu Y, Muhammad Z U H, et al. cDNA cloning, prokaryotic expression and functional analysis of 3-hydroxy-3-methylglutaryl coenzyme A reductase (HMGCR) in Pogostemon cablin[J]. Protein Expression and Purification, 2019, 163: 105454.

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1. Liu Y, Lou G, Wu W, et al. Involvement of the NF- κ B pathway in multidrug resistance induced by HBx in a hepatoma cell line[J]. Journal of viral hepatitis, 2011, 18(10): e439-e446.

2. Liu Y, Lou G, Wu W, et al. Interferon- α sensitizes HBx-expressing hepatocarcinoma cells to chemotherapeutic drugs through inhibition of HBx-mediated NF- κ B activation[J]. Virology journal, 2013, 10(1): 1-11.